# CE - DAS<sup>®</sup> UROGENICS<sup>®</sup> S – Version REF. PH00971301 CATETERE-ELETTRODO PER EMDA<sup>®</sup>

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## AVVERTENZE:

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PHYSION Srl via Fogazzaro 4/b - 41037 Mirandola (MO) - ITALY

# GB CE - DAS<sup>®</sup> UROGENICS<sup>®</sup> S – Version REF. PH00971301 CATHETER-ELECTRODE FOR EMDA<sup>®</sup>

INSTRUCTIONS FOR USE Rev. 00 of 21/11/2013

## WARNINGS:

SILICONE MEDICAL DEVICE FOR INTRAVESICAL INSTILLATIONS, WITH THE EXCEPTION OF THE BLADDER NECK, FOR USE IN A CLINICAL OR HOSPITAL ENVIRONMENT BY DOCTORS OR PARAMEDICS TRAINED TO USE THE EMDA<sup>®</sup> (ELECTRO MOTIVE DRUG ADMINISTRATION) TECHNIQUE.

DO NOT USE THE DEVICE IF THE PACKAGE IS OPENED OR DAMAGED.

THIS IS A SINGLE-USE DEVICE WHICH MUST NOT BE RISTERILISED, RECONDITIONED OR REUSED: if reused it could cause harm to the patient due to cross contaminations and/or malfunction of the device.

THIS DEVICE MUST ONLY BE USED IN COMBINATION WITH Physion's "PHYSIONIZER 30 N2", SPECIFIC AND LOW INTENSITY PULSATING POWER GENERATOR (MAX 30 MA) (REF PH03000000N) AND THE Physion dispersive electrodes (REF PG4264L400 or REF PH9473W001, with or without gel, respectively), and closely following the attached instructions.

THIS DEVICE MUST NOT BE SIMULTANEOUSLY USED WITH:

- other devices operating with electric signals, such as electrocardiographs, or which use electric current for stimulation, etc.
- other radio-frequency operated surgical devices (radio scalpels)
- other short-wave or microwave treatment devices that are running.

**STORAGE:** protected from the elements and from direct sunlight; T (-10 to +35)°C; R.H. (10 - 95)%.

**DISPOSAL:** after use, the device must be disposed of as hospital waste.

## INSTRUCTIONS:

The catheterisation procedure must absolutely be performed under ASEPTIC CONDITIONS.

1. Make sure that the spiral electrode is correctly positioned on the catheter, at the upper end.
2. Check the condition of the balloon catheter and that it inflates and deflates correctly (up to 5 ml).
3. Instill the gel in the urethra and lubricate the catheter in the conventional manner: a gel with a local anaesthetic is preferable.
4. CATHETERISE THE PATIENT laying in the supine position with the patient's legs spread slightly apart, on an appropriate surface (flat, rigid, with the pelvis not sunk into the bed) that is well lit. Inflate the balloon catheter (up to 5 ml).
5. Drain and rinse the bladder (with sterile water), then complete the drainage.
6. Instill the prescribed drug solution in the bladder.
7. Connect the red wire of the power generator to the CE-DAS<sup>®</sup> UROGENICS<sup>®</sup> S - Version catheter-electrode.
8. Connect the black wire of the power generator to the Y-cable and then to the two dispersive electrodes. Position the electrodes on the abdomen and on an area of unblemished skin.
9. Begin the therapy by exactly following the Instructions given in the User Manual included with the *PHYSIONIZER 30 N2* power generator.

**ABSOLUTE CONTRAINDICATIONS:** Do not use on patients wearing pacemakers or another active implanted medical device. Known allergies to drugs to be used. Massive intravesical bleeding denoting rupture of a large vein. Patients taking medication classified as monoamine oxidase inhibitors.

**RELATIVE CONTRAINDICATIONS:** Alcoholism - Overt psychosis - Medial lobe protruding in the bladder - Strictures of the urethra or bladder neck - Active genital-urinary tract infection (generalised inflammation).

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